

## Safety data sheet according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

- **1.1 Product identifier**
- **Trade name: LabAssay™ NEFA / Solvent B**
- **Article number:** 633-52001\_4
- **Registration number**  
Registration numbers of the individual components: see section 3.2, if applicable.
- **1.2 Relevant identified uses of the substance or mixture and uses advised against**  
No further relevant information available.
- **Application of the substance / the preparation:** Laboratory chemicals
- **Uses advised against**  
Only for professional use.  
Do not use for private purposes (household).
- **1.3 Details of the supplier of the safety data sheet**
- **Manufacturer/Supplier:**  
FUJIFILM Wako Chemicals Europe GmbH  
Fuggerstr. 12  
41468 Neuss  
Deutschland  
Tel: +49 (0) 2131-311-0  
Fax: +49 (0) 2131-311 100
- **Further information obtainable from:** msds\_eu@fujifilm.com
- **1.4 Emergency telephone number:** +49 (0) 2131-311-0 (8h-16h)

### SECTION 2: Hazards identification

- **2.1 Classification of the substance or mixture**
- **Classification according to Regulation (EC) No 1272/2008**  
The product is not classified, according to the GB CLP regulation.

---

- **2.2 Label elements**
- **Labelling according to Regulation (EC) No 1272/2008** Not applicable
- **Hazard pictograms** Not applicable
- **Signal word** Not applicable
- **Hazard statements** Not applicable
- **Additional information:**  
EUH208 Contains Mixture of 5-Chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.  
EUH210 Safety data sheet available on request.
- **2.3 Other hazards**
- **Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.

### SECTION 3: Composition/information on ingredients

- **3.2 Mixtures**
- **Description:** Mixture of substances listed below with nonhazardous additions.

(Cont. on page 2)

EN

## Safety data sheet

### according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

Trade name: LabAssay™ NEFA / Solvent B

(Cont. from page 1)

#### · Dangerous components:

55965-84-9	Mixture of 5-Chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one ⚠ Acute Tox. 3, H301; Acute Tox. 2, H310; Acute Tox. 1, H330; ⚠ Skin Corr. 1C, H314; Eye Dam. 1, H318; ⚠ Aquatic Acute 1, H400 (M=100); Aquatic Chronic 1, H410 (M=100); ⚠ Skin Sens. 1, H317 Specific concentration limits: Skin Corr. 1C; H314: C ≥ 0.6 % Skin Irrit. 2; H315: 0.06 % ≤ C < 0.6 % Eye Dam. 1; H318: C ≥ 0.6 % Eye Irrit. 2; H319: 0.06 % ≤ C < 0.6 % Skin Sens. 1; H317: C ≥ 0.0015 %	≥0.00025-<0.0015%
------------	---	-------------------

· **Additional information:** For the wording of the listed hazard phrases refer to section 16.

### SECTION 4: First aid measures

- **4.1 Description of first aid measures**
- **General information:** Immediately remove any clothing soiled by the product.
- **After inhalation:** Supply fresh air; consult doctor in case of complaints.
- **After skin contact:** Immediately wash with water and soap and rinse thoroughly.
- **After eye contact:**  
Rinse opened eye for several minutes under running water. Then consult a doctor.
- **After swallowing:** Drink plenty of water and provide fresh air. Call for a doctor immediately.
- **4.2 Most important symptoms and effects, both acute and delayed**  
No further relevant information available.
- **4.3 Indication of any immediate medical attention and special treatment needed**  
No further relevant information available.

### SECTION 5: Firefighting measures

- **5.1 Extinguishing media**
- **Suitable extinguishing agents:**  
CO<sub>2</sub>, powder or water spray. Fight larger fires with water spray or alcohol resistant foam.
- **For safety reasons unsuitable extinguishing agents:** Water with full jet
- **5.2 Special hazards arising from the substance or mixture**  
No further relevant information available.
- **5.3 Advice for firefighters**
- **Protective equipment:** No special measures required.
- **Additional information**  
Collect contaminated fire fighting water separately. It must not enter the sewage system.

### SECTION 6: Accidental release measures

- **6.1 Personal precautions, protective equipment and emergency procedures**  
Remove persons from danger area.  
Wear protective equipment. Keep unprotected persons away.

(Cont. on page 3)

## Safety data sheet according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

**Trade name: LabAssay™ NEFA / Solvent B**

(Cont. from page 2)

- **6.2 Environmental precautions:**  
Dilute with plenty of water.  
Do not allow to enter sewers/ surface or ground water.
- **6.3 Methods and material for containment and cleaning up:**  
Absorb with liquid-binding material (sand, diatomite, acid binders, universal binders, sawdust).  
Dispose contaminated material as waste according to item 13.
- **6.4 Reference to other sections**  
See Section 7 for information on safe handling.  
See Section 8 for information on personal protection equipment.  
See Section 13 for disposal information.

### SECTION 7: Handling and storage

- **7.1 Precautions for safe handling** Ensure good ventilation/exhaustion at the workplace.
- **Information about fire - and explosion protection:** No special measures required.
- **7.2 Conditions for safe storage, including any incompatibilities**
- **Storage:**
- **Requirements to be met by storerooms and receptacles:**  
Store only in the original receptacle.  
Store in cool, dry conditions in well sealed receptacles.
- **Information about storage in one common storage facility:** Store away from foodstuffs.
- **Further information about storage conditions:**  
Protect from frost.  
2-10 °C
- **Storage class:** 12
- **7.3 Specific end use(s)** No further relevant information available.

### SECTION 8: Exposure controls/personal protection

- **8.1 Control parameters**
- **Ingredients with limit values that require monitoring at the workplace:**  
The product does not contain any relevant quantities of materials with critical values that have to be monitored at the workplace.
- **Additional information:** The lists valid during the making were used as basis.
- **8.2 Exposure controls**
- **Appropriate engineering controls** No further data; see item 7.
- **Individual protection measures, such as personal protective equipment**
- **General protective and hygienic measures:**  
Wash hands before breaks and at the end of work.  
Avoid contact with the eyes and skin.  
Do not eat, drink, smoke or sniff while working.  
Storing food in the working area is prohibited.
- **Respiratory protection:** Not necessary if room is well-ventilated.
- **Hand protection**



Protective gloves

(Cont. on page 4)

EN

## Safety data sheet

### according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

Trade name: LabAssay™ NEFA / Solvent B

(Cont. from page 3)

- **Material of gloves**  
Use chemical resistant gloves. In case of prolonged immersion or frequently repeated contact use gloves made of the materials: nitrile rubber (thickness  $\geq$  0.38 mm, breakthrough time > 480 min).
- **Eye/face protection**



Safety glasses

The protective safety glasses to be used must comply with the specifications of the EC directive 89/686/EEC and the resultant standard EN166.

- **Body protection:** Protective work clothing

## SECTION 9: Physical and chemical properties

### · 9.1 Information on basic physical and chemical properties

- |   |                                    |
|---|------------------------------------|
| · <b>General Information</b>                                      |                                    |
| · <b>Physical state</b>   | Fluid                              |
| · <b>Colour:</b>  | According to product specification |
| · <b>Odour:</b>   | Characteristic                     |
| · <b>Odour threshold:</b>   | Not determined.                    |
| · <b>Melting point/freezing point:</b>                            | Undetermined.                      |
| · <b>Boiling point or initial boiling point and boiling range</b> | Undetermined.                      |
| · <b>Flammability</b>   | Not applicable.                    |
| · <b>Lower and upper explosion limit</b>                          |                                    |
| · <b>Lower:</b>   | Not determined.                    |
| · <b>Upper:</b>   | Not determined.                    |
| · <b>Flash point:</b>   | Not applicable.                    |
| · <b>Ignition temperature:</b>                                    | Not determined.                    |
| · <b>Decomposition temperature:</b>                               | Not determined.                    |
| · <b>pH</b>   | Not determined.                    |
| · <b>Viscosity:</b>   |                                    |
| · <b>Kinematic viscosity</b>                                      | Not determined.                    |
| · <b>Dynamic:</b>   | Not determined.                    |
| · <b>Solubility</b>   |                                    |
| · <b>water:</b>   | Fully miscible.                    |
| · <b>Partition coefficient n-octanol/water (log value)</b>        | Not determined.                    |
| · <b>Vapour pressure:</b>   | Not determined.                    |
| · <b>Density and/or relative density</b>                          |                                    |
| · <b>Density:</b>   | Not determined.                    |
| · <b>Relative density</b>   | Not determined.                    |
| · <b>Vapour density</b>   | Not determined.                    |

### · 9.2 Other information

- |                                     |                              |
|-------------------------------------|------------------------------|
| · <b>Appearance:</b>                |                              |
| · <b>Form:</b>                      | Liquid                       |
| · <b>Auto-ignition temperature:</b> | Product is not selfigniting. |

(Cont. on page 5)

## Safety data sheet

### according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

Trade name: LabAssay™ NEFA / Solvent B

(Cont. from page 4)

- |                                |                 |
|--------------------------------|-----------------|
| · <b>Explosive properties:</b> | Not determined. |
| · <b>Evaporation rate</b>      | Not determined. |

· **Information with regard to physical hazard classes**

- |  |                |
|--|----------------|
| · <b>Explosives</b>  | Not applicable |
| · <b>Flammable gases</b>   | Not applicable |
| · <b>Aerosols</b>  | Not applicable |
| · <b>Oxidising gases</b>   | Not applicable |
| · <b>Gases under pressure</b>  | Not applicable |
| · <b>Flammable liquids</b>   | Not applicable |
| · <b>Flammable solids</b>  | Not applicable |
| · <b>Self-reactive substances and mixtures</b>                                     | Not applicable |
| · <b>Pyrophoric liquids</b>  | Not applicable |
| · <b>Pyrophoric solids</b>   | Not applicable |
| · <b>Self-heating substances and mixtures</b>                                      | Not applicable |
| · <b>Substances and mixtures, which emit flammable gases in contact with water</b> | Not applicable |
| · <b>Oxidising liquids</b>   | Not applicable |
| · <b>Oxidising solids</b>  | Not applicable |
| · <b>Organic peroxides</b>   | Not applicable |
| · <b>Corrosive to metals</b>   | Not applicable |
| · <b>Desensitised explosives</b>   | Not applicable |

### SECTION 10: Stability and reactivity

- **10.1 Reactivity** No further relevant information available.
- **10.2 Chemical stability**
- **Thermal decomposition / conditions to be avoided:**  
No decomposition if used according to specifications.
- **10.3 Possibility of hazardous reactions** No dangerous reactions known.
- **10.4 Conditions to avoid**  
See "Fire Fighting measures".  
Protect from heat and direct sunlight.
- **10.5 Incompatible materials:** No further relevant information available.
- **10.6 Hazardous decomposition products:** No dangerous decomposition products known.

### SECTION 11: Toxicological information

- **11.2 Information on other hazards**
- **Endocrine disrupting properties**

None of the ingredients is listed.

EN

(Cont. on page 6)

## Safety data sheet according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

Trade name: LabAssay™ NEFA / Solvent B

(Cont. from page 5)

### SECTION 12: Ecological information

- **12.1 Toxicity**
- **Aquatic toxicity:** No further relevant information available.
- **12.2 Persistence and degradability** No further relevant information available.
- **12.3 Bioaccumulative potential** No further relevant information available.
- **12.4 Mobility in soil** No further relevant information available.
- **12.5 Results of PBT and vPvB assessment**
- **PBT:** Not applicable.
- **vPvB:** Not applicable.
- **12.6 Endocrine disrupting properties**  
The product does not contain substances with endocrine disrupting properties.
- **12.7 Other adverse effects**
- **Additional ecological information:**
- **General notes:**  
Do not allow undiluted product or large quantities of it to reach ground water, water course or sewage system.

### SECTION 13: Disposal considerations

- **13.1 Waste treatment methods**
- **Recommendation**  
Must not be disposed together with household garbage. Do not allow product to reach sewage system.  
Contact waste processors for recycling information.
- **Uncleaned packaging:**
- **Recommendation:**  
Disposal must be made according to official regulations.  
Packagings that may not be cleansed are to be disposed of in the same manner as the product.
- **Recommended cleansing agents:** Water, if necessary together with cleansing agents.

### SECTION 14: Transport information

- |   |                 |
|---|-----------------|
| · <b>14.1 UN number or ID number</b><br>· <b>ADR, IMDG, IATA</b>                            | Not applicable  |
| · <b>14.2 UN proper shipping name</b><br>· <b>ADR, IMDG, IATA</b>                           | Not applicable  |
| · <b>14.3 Transport hazard class(es)</b><br>· <b>ADR, ADN, IMDG, IATA</b><br>· <b>Class</b> | Not applicable  |
| · <b>14.4 Packing group</b><br>· <b>ADR, IMDG, IATA</b>                                     | Not applicable  |
| · <b>14.5 Environmental hazards:</b>  | Not applicable. |

(Cont. on page 7)

## Safety data sheet according to 1907/2006/EC, Article 31

Printing date 09.01.2023

Version number 1.0

Revision: 10.05.2022

Trade name: LabAssay™ NEFA / Solvent B

(Cont. from page 6)

- |   |                 |
|---|-----------------|
| · <b>14.6 Special precautions for user</b>                            | Not applicable. |
| · <b>14.7 Maritime transport in bulk according to IMO instruments</b> | Not applicable. |
| · <b>UN "Model Regulation":</b>                                       | Not applicable  |

### SECTION 15: Regulatory information

- **15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture**
- **Directive 2012/18/EU**
- **Named dangerous substances - ANNEX I** None of the ingredients is listed.
- **15.2 Chemical safety assessment:** A Chemical Safety Assessment has not been carried out.

### SECTION 16: Other information

This information is based on our present knowledge. However, this shall not constitute a guarantee for any specific product features and shall not establish a legally valid contractual relationship.

- **Relevant phrases**

- H301 Toxic if swallowed.
- H310 Fatal in contact with skin.
- H314 Causes severe skin burns and eye damage.
- H317 May cause an allergic skin reaction.
- H318 Causes serious eye damage.
- H330 Fatal if inhaled.
- H400 Very toxic to aquatic life.
- H410 Very toxic to aquatic life with long lasting effects.

- **Creation Date** 09.05.2022

- **Abbreviations and acronyms:**

- ADR: Accord relatif au transport international des marchandises dangereuses par route (European Agreement Concerning the International Carriage of Dangerous Goods by Road)
- IMDG: International Maritime Code for Dangerous Goods
- IATA: International Air Transport Association
- GHS: Globally Harmonised System of Classification and Labelling of Chemicals
- EINECS: European Inventory of Existing Commercial Chemical Substances
- ELINCS: European List of Notified Chemical Substances
- CAS: Chemical Abstracts Service (division of the American Chemical Society)
- PBT: Persistent, Bioaccumulative and Toxic
- vPvB: very Persistent and very Bioaccumulative
- Acute Tox. 3: Acute toxicity – Category 3
- Acute Tox. 2: Acute toxicity – Category 2
- Acute Tox. 1: Acute toxicity – Category 1
- Skin Corr. 1C: Skin corrosion/irritation – Category 1C
- Eye Dam. 1: Serious eye damage/eye irritation – Category 1
- Skin Sens. 1: Skin sensitisation – Category 1
- Aquatic Acute 1: Hazardous to the aquatic environment - acute aquatic hazard – Category 1
- Aquatic Chronic 1: Hazardous to the aquatic environment - long-term aquatic hazard – Category 1

- **\* Data modified compared to previous version.**