

SAFETY DATA SHEET

According to JIS Z 7253:2019
Revision Date 1-Jul-2023
Version 2

Section 1: PRODUCT AND COMPANY IDENTIFICATION

Product name	LBIS Rat LH ELISA Kit (S-type)
Product code	636-23921

Manufacturer	FUJIFILM Wako Pure Chemical Corporation 1-2 Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Phone: +81-6-6203-3741 Facsimile: +81-6-6203-2029
Supplier	FUJIFILM Wako Pure Chemical Corporation 1-2 Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Phone: +81-6-6203-3741 Facsimile: +81-6-6203-2029
Emergency telephone number	+81-6-6203-3741 / +81-3-3270-8571
Recommended uses and restrictions on use	For research use only

Section 2: HAZARDS IDENTIFICATION

GHS classification

Classification of the substance or mixture

Acute toxicity - Inhalation (Vapors)	Category 4
Skin corrosion/irritation	Category 1
Serious eye damage/eye irritation	Category 1
Skin sensitization	Category 1
Specific target organ toxicity (single exposure)	Category 1
Category 1 respiratory system	
Specific target organ toxicity (repeated exposure)	Category 1
Category 1 respiratory system	

Pictograms



Signal word

Danger

Hazard statements

- H314 - Causes severe skin burns and eye damage
- H318 - Causes serious eye damage
- H332 - Harmful if inhaled
- H317 - May cause an allergic skin reaction
- H370 - Causes damage to the following organs: respiratory system
- H372 - Causes damage to the following organs through prolonged or repeated exposure: respiratory system

Precautionary statements-(Prevention)

- Use only outdoors or in a well-ventilated area
- Do not breathe dust/ fume/ gas/ mist/ vapors/spray
- Wash face, hands and any exposed skin thoroughly after handling
- Wear protective gloves/protective clothing/eye protection/face protection
- Contaminated work clothing should not be allowed out of the workplace
- Do not eat, drink or smoke when using this product

Precautionary statements-(Response)

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- Immediately call a POISON CENTER or doctor/physician
- IF ON SKIN (or hair): Remove/ Take off immediately all contaminated clothing. Rinse skin with water/shower.
- Wash contaminated clothing before reuse.
- If skin irritation or rash occurs: Get medical advice/ attention
- IF INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
- Call a POISON CENTER or doctor /physician if you feel unwell.
- IF SWALLOWED: Rinse mouth. DO NOT induce vomiting

Precautionary statements-(Storage)

- Store locked up.

Precautionary statements-(Disposal)

- Dispose of contents/container to an approved waste disposal plant

Others

Other hazards Not available

Section 3: COMPOSITION/INFORMATION ON INGREDIENTS**Single Substance or Mixture**

Kit (Set of mixtures)

Chemical Name	Weight-%	Molecular weight	ENCS	ISHL No.	CAS RN
Antibody-coated 96 well-plate	-	N/A	N/A	N/A	N/A
Rat standard LH	-	N/A	N/A	N/A	N/A
Assay buffer	-	N/A	N/A	N/A	N/A
Biotinylated anti-LH α antibody	-	N/A	N/A	N/A	N/A
HRP-conjugated streptavidin	-	N/A	N/A	N/A	N/A
Chromogen (TMB): 3,3',5,5' - Tetramethyl-benzidine	-	N/A	N/A	N/A	N/A
Sample dilution buffer	-	N/A	N/A	N/A	N/A
Stop solution: Sulfuric Acid	-	N/A	N/A	N/A	N/A
Wash stock solution (10X): Phosphate Buffer	-	N/A	N/A	N/A	N/A

Impurities and/or Additives : Not applicable

Substances Remarks: Sulfuric Acid 6.9 %, Poly (oxyethylene) sorbitan monolaurate <1 %, 2-Methyl-2H-isothiazol-3-one <0.2 %

Section 4: FIRST AID MEASURES**Inhalation**

Remove to fresh air. If symptoms persist, call a physician.

Skin contact

Wash off immediately with soap and plenty of water. If symptoms persist, call a physician.

Eye contact

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy

to do. Continue rinsing. Immediate medical attention is required.

Ingestion

Rinse mouth. Never give anything by mouth to an unconscious person. Call a physician or poison control center immediately. Do not induce vomiting without medical advice.

Protection of first-aiders

Use personal protective equipment as required.

Section 5: FIRE FIGHTING MEASURES

Suitable extinguishing media

Water spray (fog), Carbon dioxide (CO₂), Foam, Extinguishing powder, Sand

Unsuitable extinguishing media

No information available

Specific hazards arising from the chemical product

Thermal decomposition can lead to release of irritating and toxic gases and vapors.

Special extinguishing method

No information available

Special protective actions for fire-fighters

Use personal protective equipment as required. Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear.

Section 6: ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

For indoor, provide adequate ventilation process until the end of working. Deny unnecessary entry other than the people involved by, for example, using a rope. While working, wear appropriate protective equipments to avoid adhering it on skin, or inhaling the gas. Work from windward, and retract the people downwind.

Environmental precautions

To be careful not discharged to the environment without being properly handled waste water contaminated.

Methods and materials for contaminant and methods and materials for cleaning up

Absorb dry sand, earth, sawdust and the waste. Collect empty container that can be sealed.

Recovery, neutralization

No information available

Secondary disaster prevention measures

Clean contaminated objects and areas thoroughly observing environmental regulations.

Section 7: HANDLING AND STORAGE

Handling

Technical measures

Avoid contact with oxidizing and reducing agents. Avoid contact with alkaline substances. Use with local exhaust ventilation.

Precautions

Do not rough handling containers, such as upsetting, falling, giving a shock, and dragging. Prevent leakage, overflow, and scattering. Not to generate steam and dust in vain. Seal the container after use. After handling, wash hands and face, and then gargle. In places other than those specified, should not be smoking or eating and drinking. Should not be brought contaminated protective equipment and gloves to rest stops. Deny unnecessary entry of non-emergency personnel to the handling area.

Safety handling precautions

Use personal protective equipment as required. Avoid contact with skin, eyes or clothing.

Storage

Safe storage conditions

Storage conditions Store away from sunlight in a cool (2 °C -8 °C) well-ventilated dry place.

Safe packaging material No information available

Incompatible substances Strong oxidizing agents, Reducing agent, Alkali

Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

Engineering controls

In case of indoor workplace, seal the source or use a local exhaust system. Provide the safety shower facility, and hand- and eye-wash facility. And display their position clearly.

Exposure limits

Chemical Name	JSOH (Japan)	ISHL (Japan)	ACGIH
Sulfuric Acid 7664-93-9	1 mg/m ³	N/A	TWA 0.2 mg/m ³

Personal protective equipment

Respiratory protection Gas mask for acidic gas
Hand protection Impermeable protective gloves
Eye protection protective eyeglasses or chemical safety goggles
Skin and body protection Long-sleeved work clothes

General hygiene considerations

Handle in accordance with good industrial hygiene and safety practice.

Section 9: PHYSICAL AND CHEMICAL PROPERTIES

Form

Appearance	Kit (Set of mixtures)
Odor	No data available
Melting point/freezing point	No data available
Boiling point, initial boiling point and boiling range	No data available
Flammability	No data available
Evaporation rate:	No data available
Flammability (solid, gas):	No data available
Upper/lower flammability or explosive limits	
Upper :	No data available
Lower :	No data available
Flash point	No data available
Auto-ignition temperature:	No data available
Decomposition temperature:	No data available
pH	No data available
Viscosity (coefficient of viscosity)	No data available
Dynamic viscosity	No data available
Solubilities	No data available
n-Octanol/water partition coefficient: (log Pow)	No data available
Vapor pressure	No data available
Specific Gravity/ Relative density	No data available
Vapor density	No data available
Particle characteristics	No data available

Section 10: STABILITY AND REACTIVITY

Stability

Reactivity No data available
Chemical stability Stable under recommended storage conditions.

Hazardous reactions

None under normal processing

Conditions to avoid

Extremes of temperature and direct sunlight

Incompatible materials

Strong oxidizing agents, Reducing agent, Alkali

Hazardous decomposition productsCarbon monoxide (CO), Carbon dioxide (CO₂), Nitrogen oxides (NO_x), Sulfur oxides (SO_x), Phosphorus oxide**Section 11: TOXICOLOGICAL INFORMATION****Acute toxicity**

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Sulfuric Acid	2140 mg/kg (Rat)	N/A	347 ppm (Rat) 4 h
Poly (oxyethylene) sorbitan monolaurate	37000 mg/kg (Rat) 36700 µL/kg (Rat)	N/A	N/A
2-Methyl-2H-isothiazol-3-one	120 mg/kg (Rat)	200 mg/kg (Rabbit)	0.11 mg/L (Rat) 4 h

Chemical Name	Acute toxicity -oral-source information	Acute toxicity -dermal-source information	Acute toxicity -inhalation gas- source information
Sulfuric Acid	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.

Chemical Name	Acute toxicity -inhalation vapor- source information	Acute toxicity -inhalation dust- source information	Acute toxicity -inhalation mist- source information
Sulfuric Acid	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.	Based on the NITE GHS classification results.

Skin irritation/corrosion

Chemical Name	Skin corrosion/irritation source information
Sulfuric Acid	Based on the NITE GHS classification results.

Serious eye damage/ irritation

Chemical Name	Serious eye damage/irritation source information
Sulfuric Acid	Based on the NITE GHS classification results.

Respiratory or skin sensitization

Chemical Name	Respiratory or Skin sensitization source information
Sulfuric Acid	Based on the NITE GHS classification results.

Reproductive cell mutagenicity

Chemical Name	germ cell mutagenicity source information
Sulfuric Acid	Based on the NITE GHS classification results.

Carcinogenicity

Chemical Name	Carcinogenicity source information
Sulfuric Acid	Based on the NITE GHS classification results.

Chemical Name	NTP	IARC	ACGIH	JSOH (Japan)
Sulfuric Acid 7664-93-9	-	Group 1	A2	-

Reproductive toxicity

Chemical Name	Reproductive toxicity source information
Sulfuric Acid	Based on the NITE GHS classification results.

STOT-single exposure

Chemical Name	STOT -single exposure- source information
---------------	---

Sulfuric Acid	Based on the NITE GHS classification results.
---------------	---

STOT-repeated exposure

Chemical Name	STOT -repeated exposure- source information
Sulfuric Acid	Based on the NITE GHS classification results.

Aspiration hazard

Chemical Name	Aspiration Hazard source information
Sulfuric Acid	Based on the NITE GHS classification results.

Section 12: ECOLOGICAL INFORMATION

Ecotoxicity

Chemical Name	Algae/aquatic plants	Fish	Crustacea
Sulfuric Acid	N/A	LC50: <i>Lepomis macrochirus</i> 16 - 28mg/L 96 h	LC50: <i>Daphnia magna</i> 29 mg/L 24 h
2-Methyl-2H-isothiazol-3-one	N/A	LC50: <i>Oncorhynchus mykiss</i> 0.07 mg/L 96h	EC50: <i>Daphnia magna</i> 0.18 mg/L 48h

Other data

Chemical Name	Short-term (acute) hazardous to the aquatic environment source information	Long-term (chronic) hazardous to the aquatic environment source information
Sulfuric Acid	Based on the NITE GHS classification results	Based on the NITE GHS classification results

Persistence and degradability	No information available
Bioaccumulative potential	No information available
Mobility in soil	No information available
Hazard to the ozone layer	No information available

Section 13: DISPOSAL CONSIDERATIONS

Waste from residues

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Contaminated container and contaminated packaging

Disposal should be in accordance with applicable regional, national and local laws and regulations.

Section 14: TRANSPORT INFORMATION

ADR/RID

UN number	UN3264
Proper shipping name:	Corrosive liquid, acidic, inorganic, n.o.s. (Diluted Sulfuric Acid)
UN classification	8
Subsidiary hazard class	
Packing group	II
Marine pollutant	Not applicable

IMDG

UN number	UN3264
Proper shipping name:	Corrosive liquid, acidic, inorganic, n.o.s. (Diluted Sulfuric Acid)
UN classification	8
Subsidiary hazard class	
Packing group	II
Marine pollutant (Sea)	Not applicable

Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code No information available

IATA

UN number UN3264
 Proper shipping name: Corrosive liquid, acidic, inorganic, n.o.s. (Diluted Sulfuric Acid)
 UN classification 8
 Subsidiary hazard class
 Packing group II
 Environmentally Hazardous Substance Not applicable

Section 15: REGULATORY INFORMATION

International Inventories

EINECS/ELINCS -
 TSCA -

Japanese regulations

Fire Service Act Not applicable
Poisonous and Deleterious Substances Control Law Not applicable
Industrial Safety and Health Act Harmful Substances Whose Names Are to be Indicated on the Label (Law Art.57, Para.1, Enforcement Order Art.18) Group 3 Specified Chemical Substance, (Ordinance on Prevention of Hazards Due to Specified Chemical Substances Art.2 Para.1, Item 6) Notifiable Substances (Law Art.57-2, Enforcement Order Art.18-2 Attached Table No.9)No.613 Priority Assessment Chemical Substances (Law Article 2, Para.5)
Act on the Evaluation of Chemical Substances and Regulation of Their Manufacture, etc
Regulations for the carriage and storage of dangerous goods in ship Corrosive Substances (Ordinance Art.3, Ministry of Transportation Ordinance Regarding Transport by Ship and Storage, Attached Table 1)
Civil Aeronautics Law Corrosive Substances (Ordinance Art.194, MITL Notification for Air Transportation of Explosives etc., Attached Table 1)
Marine Pollution Prevention Law Enforcement ordinance Appendix No.1 Noxious liquid substance category Y
Pollutant Release and Transfer Register Law Not applicable
Water Pollution Control Act Specified substances (Law Art.2 Para.4, Enforcement Order Art.3-3)
Air Pollution Control Law Specified Substances

Industrial Safety and Health Law

Law Name	Chemical Name in Regulation	Ordinance Number	Weight %
Notifiable Substances (Law Art.57-2, Enforcement Order Art.18-2 Attached Table No.9, and Law Art.56-1)	Sulfuric acid	613	6.9

Section 16: OTHER INFORMATION

Key literature references and sources NITE: National Institute of Technology and Evaluation (JAPAN)

for data etc.

<http://www.safe.nite.go.jp/japan/db.html>
IATA dangerous Goods Regulations
RTECS: Registry of Toxic Effects of Chemical Substances
Japan Industrial Safety and Health Association GHS Model SDS
Dictionary of Synthetic Organic Chemistry, SSOCJ, Koudansha
Scientific Co.Ltd.
Chemical Dictionary, Kyouritsu Publishing Co., Ltd.
etc

Disclaimer

This SDS is according to JIS Z 7253: 2019. The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text. GHS Classification is according to JIS Z7252(2019). *JIS: Japanese Industrial Standards

End of Safety Data Sheet