

## SAFETY DATA SHEET

According to JIS Z 7253:2019  
**Revision Date** 1-Jul-2023  
 Version 3

## Section 1: PRODUCT AND COMPANY IDENTIFICATION

<b>Product name</b>	LBIS Mouse/Rat Insulin ELISA Kit (Luminescent type)
<b>Product code</b>	637-54721

<b>Manufacturer</b>	FUJIFILM Wako Pure Chemical Corporation 1-2 Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Phone: +81-6-6203-3741 Facsimile: +81-6-6203-2029
<b>Supplier</b>	FUJIFILM Wako Pure Chemical Corporation 1-2 Doshomachi 3-Chome, Chuo-ku, Osaka 540-8605, Japan Phone: +81-6-6203-3741 Facsimile: +81-6-6203-2029
<b>Emergency telephone number</b>	+81-6-6203-3741 / +81-3-3270-8571
<b>Recommended uses and restrictions on use</b>	For research use only

## Section 2: HAZARDS IDENTIFICATION

**GHS classification****Classification of the substance or mixture**

**Serious eye damage/eye irritation**

Category 2B

**Skin sensitization**

Category 1

**Pictograms****Signal word**

Warning

**Hazard statements**

H320 - Causes eye irritation

H317 - May cause an allergic skin reaction

**Precautionary statements-(Prevention)**

- Wash face, hands and any exposed skin thoroughly after handling
- Avoid breathing dust/ fume/ gas/ mist/ vapors/ spray
- Contaminated work clothing should not be allowed out of the workplace
- Protective gloves

**Precautionary statements-(Response)**

- IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
- If eye irritation persists: Get medical advice/ attention.
- IF ON SKIN: Wash with plenty of soap and water
- If skin irritation or rash occurs: Get medical advice/ attention
- Wash contaminated clothing before reuse.

**Precautionary statements-(Storage)**

- Not applicable

**Precautionary statements-(Disposal)**

- Dispose of contents/container to an approved waste disposal plant

**Others****Other hazards**

Not available

**Section 3: COMPOSITION/INFORMATION ON INGREDIENTS****Single Substance or Mixture**

Kit (Set of mixtures)

Chemical Name	Weight-%	Molecular weight	ENCS	ISHL No.	CAS RN
Anti-Insulin antibody-coated plate	-	N/A	N/A	N/A	N/A
Standard mouse/rat Insulin solution	-	N/A	N/A	N/A	N/A
Buffer solution	-	N/A	N/A	N/A	N/A
HRP-conjugated anti-insulin antibody	-	N/A	N/A	N/A	N/A
Luminescent Reagent 1	-	N/A	N/A	N/A	N/A
Luminescent Reagent 2	-	N/A	N/A	N/A	N/A
Wash stock solution (10X)	-	N/A	N/A	N/A	N/A

**Impurities and/or Additives :**

Not applicable

**Hazardous Component**Polyoxyethylene (20) Sorbitan Monolaurate 0.5 %,  
2-Methyl-2H-isothiazol-3-one <0.2 %**Section 4: FIRST AID MEASURES****Inhalation**

Remove to fresh air. If symptoms persist, call a physician.

**Skin contact**

Wash off immediately with soap and plenty of water. If symptoms persist, call a physician.

**Eye contact**

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediate medical attention is required.

**Ingestion**

Rinse mouth. Never give anything by mouth to an unconscious person. Call a physician or poison control center immediately. Do not induce vomiting without medical advice.

**Protection of first-aiders**

Use personal protective equipment as required.

**Section 5: FIRE FIGHTING MEASURES****Suitable extinguishing media**

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment

**Unsuitable extinguishing media**

No information available

**Specific hazards arising from the chemical product**

Thermal decomposition can lead to release of irritating and toxic gases and vapors.

**Special extinguishing method**

No information available

**Special protective actions for fire-fighters**

Use personal protective equipment as required. Firefighters should wear self-contained breathing apparatus and full firefighting turnout gear.

## Section 6: ACCIDENTAL RELEASE MEASURES

### Personal precautions, protective equipment and emergency procedures

For indoor, provide adequate ventilation process until the end of working. Deny unnecessary entry other than the people involved by, for example, using a rope. While working, wear appropriate protective equipments to avoid adhering it on skin, or inhaling the gas. Work from windward, and retract the people downwind.

### Environmental precautions

To be careful not discharged to the environment without being properly handled waste water contaminated.

### Methods and materials for contaminant and methods and materials for cleaning up

Absorb dry sand, earth, sawdust and the waste. Collect empty container that can be sealed.

### Recovery, neutralization

No information available

### Secondary disaster prevention measures

Clean contaminated objects and areas thoroughly observing environmental regulations.

## Section 7: HANDLING AND STORAGE

### Handling

#### Technical measures

Avoid contact with strong oxidizing agents. Use with local exhaust ventilation.

#### Precautions

Do not rough handling containers, such as upsetting, falling, giving a shock, and dragging. Prevent leakage, overflow, and scattering. Not to generate steam and dust in vain. Seal the container after use. After handling, wash hands and face, and then gargle. In places other than those specified, should not be smoking or eating and drinking. Should not be brought contaminated protective equipment and gloves to rest stops. Deny unnecessary entry of non-emergency personnel to the handling area.

#### Safety handling precautions

Use personal protective equipment as required. Avoid contact with skin, eyes or clothing.

#### Storage

##### Safe storage conditions

**Storage conditions** Store away from sunlight in a cool (2 °C -8 °C) well-ventilated dry place.

**Safe packaging material** No information available

**Incompatible substances** Strong oxidizing agents

## Section 8: EXPOSURE CONTROLS/PERSONAL PROTECTION

### Engineering controls

In case of indoor workplace, seal the source or use a local exhaust system. Provide the safety shower facility, and hand- and eye-wash facility. And display their position clearly.

### Exposure limits

This product, as supplied, does not contain any hazardous materials with occupational exposure limits established by the region specific regulatory bodies.

### Personal protective equipment

**Respiratory protection** Protective mask

**Hand protection** Protection gloves

**Eye protection** protective eyeglasses or chemical safety goggles

**Skin and body protection** Long-sleeved work clothes

### General hygiene considerations

Handle in accordance with good industrial hygiene and safety practice.

## Section 9: PHYSICAL AND CHEMICAL PROPERTIES

<b>Form</b>	
<b>Appearance</b>	Kit (Set of mixtures)
<b>Odor</b>	No data available
<b>Melting point/freezing point</b>	No data available
<b>Boiling point, initial boiling point and boiling range</b>	No data available
<b>Flammability</b>	No data available
<b>Evaporation rate:</b>	No data available
<b>Flammability (solid, gas):</b>	No data available
<b>Upper/lower flammability or explosive limits</b>	
<b>Upper :</b>	No data available
<b>Lower :</b>	No data available
<b>Flash point</b>	No data available
<b>Auto-ignition temperature:</b>	No data available
<b>Decomposition temperature:</b>	No data available
<b>pH</b>	No data available
<b>Viscosity (coefficient of viscosity)</b>	No data available
<b>Dynamic viscosity</b>	No data available
<b>Solubilities</b>	No data available
<b>n-Octanol/water partition coefficient: (log Pow)</b>	No data available
<b>Vapor pressure</b>	No data available
<b>Specific Gravity / Relative density</b>	No data available
<b>Vapor density</b>	No data available
<b>Particle characteristics</b>	No data available

## Section 10: STABILITY AND REACTIVITY

<b>Stability</b>	
<b>Reactivity</b>	No data available
<b>Chemical stability</b>	Stable under recommended storage conditions.
<b>Hazardous reactions</b>	
None under normal processing	
<b>Conditions to avoid</b>	
Extremes of temperature and direct sunlight	
<b>Incompatible materials</b>	
Strong oxidizing agents	
<b>Hazardous decomposition products</b>	
Carbon monoxide (CO), Carbon dioxide (CO <sub>2</sub> ), Nitrogen oxides (NO <sub>x</sub> ), Phosphorus oxide	

## Section 11: TOXICOLOGICAL INFORMATION

### Acute toxicity

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Poly (oxyethylene) sorbitan monolaurate	37000 mg/kg (Rat) 36700 µL/kg (Rat)	N/A	N/A
2-Methyl-2H-isothiazol-3-one	120 mg/kg (Rat)	200 mg/kg (Rabbit)	0.11 mg/L (Rat) 4 h

<b>Skin irritation/corrosion</b>	No data available
<b>Serious eye damage/ irritation</b>	No data available
<b>Respiratory or skin sensitization</b>	No data available
<b>Reproductive cell mutagenicity</b>	No data available

<b>Carcinogenicity</b>	No data available
<b>Reproductive toxicity</b>	No data available
<b>STOT-single exposure</b>	No data available
<b>STOT-repeated exposure</b>	No data available
<b>Aspiration hazard</b>	No data available

### Section 12: ECOLOGICAL INFORMATION

#### Ecotoxicity

Chemical Name	Algae/aquatic plants	Fish	Crustacea
2-Methyl-2H-isothiazol-3-one	N/A	LC50: <i>Oncorhynchus mykiss</i> 0.07 mg/L, 96h	EC50: <i>Daphnia magna</i> 0.18 mg/L, 48h

**Other data** No data available

<b>Persistence and degradability</b>	No information available
<b>Bioaccumulative potential</b>	No information available
<b>Mobility in soil</b>	No information available
<b>Hazard to the ozone layer</b>	No information available
<b>Mobility</b>	

### Section 13: DISPOSAL CONSIDERATIONS

#### Waste from residues

Disposal should be in accordance with applicable regional, national and local laws and regulations.

#### Contaminated container and contaminated packaging

Disposal should be in accordance with applicable regional, national and local laws and regulations.

### Section 14: TRANSPORT INFORMATION

<b>ADR/RID</b>	Not regulated
<b>UN number</b>	-
<b>Proper shipping name:</b>	
<b>UN classification</b>	
<b>Subsidiary hazard class</b>	
<b>Packing group</b>	
<b>Marine pollutant</b>	Not applicable
<b>IMDG</b>	Not regulated
<b>UN number</b>	-
<b>Proper shipping name:</b>	
<b>UN classification</b>	
<b>Subsidiary hazard class</b>	
<b>Packing group</b>	
<b>Marine pollutant (Sea)</b>	Not applicable
<b>Transport in bulk according to Annex II of MARPOL 73/78 and the IBC Code</b>	No information available
<b>IATA</b>	Not regulated
<b>UN number</b>	-
<b>Proper shipping name:</b>	
<b>UN classification</b>	
<b>Subsidiary hazard class</b>	

Packing group  
Environmentally Hazardous Substance Not applicable

### Section 15: REGULATORY INFORMATION

#### International Inventories

EINECS/ELINCS -  
TSCA -

#### Japanese regulations

Fire Service Act	Not applicable
Poisonous and Deleterious Substances Control Law	Not applicable
Industrial Safety and Health Act	Not applicable
Act on the Evaluation of Chemical Substances and Regulation of Their Manufacture, etc	Priority Assessment Chemical Substances (Law Article 2, Para.5)
Regulations for the carriage and storage of dangerous goods in ship	Not applicable
Civil Aeronautics Law	Not applicable
Pollutant Release and Transfer Register Law	Not applicable

### Section 16: OTHER INFORMATION

**Key literature references and sources for data etc.** NITE: National Institute of Technology and Evaluation (JAPAN)  
<http://www.safe.nite.go.jp/japan/db.html>  
IATA dangerous Goods Regulations  
RTECS: Registry of Toxic Effects of Chemical Substances  
Japan Industrial Safety and Health Association GHS Model SDS  
Dictionary of Synthetic Organic Chemistry, SSOCJ, Koudansha Scientific Co.Ltd.  
Chemical Dictionary, Kyouritsu Publishing Co., Ltd.  
etc

#### **Disclaimer**

This SDS is according to JIS Z 7253: 2019. The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text. GHS Classification is according to JIS Z7252(2019). \*JIS: Japanese Industrial Standards

**End of Safety Data Sheet**